



**CL Probe CC1N Electrode Cable Connections to
AGSELECT B2 - ESELECT B2 - ESELECT M2**

ACV0100404 5 PIN CC1N TO ESELECT M2



ACV0000404 4 PIN CC1N to ESELECT M2



Previous Generation Controllers:

AG Select & eSelect B2 with old style 4 Pin CC1N Probe Connector – Original Cable Code SCV0001003 with BNC and Black Elbow Connections. This cable can be modified by removal of the BNC and Black Elbow, then trim back to utilise the White, Brown, Green & Yellow wires.

AG Select & eSelect B2 with new 5 Pin CC1N Probe Connector - Use SCV0101003 5 PIN to AG SELECT or ESELECT B2 to suit new type 5 pin electrode.